ANSWERS TO EXERCISE IN SESSION 3.

1. MAKE OBSERVATIONS. All living things have DNA, it is inside the cells.
2. COME UP WITH A QUESTION BASED ON OBSERVATIONS: Can we extract DNA from a strawberry? How can we do it?
3. MAKE A HYPOTHESIS: We think that we can extract DNA from a strawberry because it comes from a plant, and a plant is a living thing.
4. CONDUCT EXPERIMENTS:
 - a. Smash the strawberries inside a plastic zip lock bag.
 - b. Make the DNA extraction liquid (lysis): 2 tea spoons of detergent, 1 tea spoon of salt and half a cup of water.
 - c. Pour the liquid inside the bag with the strawberries and mix it.
 - d. Filter the mixture thanks to a coffee filter.
 - e. Pour the same amount of alcohol and wait.
 - f. A cloudy substance will appear on top. It is DNA.
5. RECORD YOUR RESULTS AND DRAW A CONCLUSION. We learned that we can extract DNA from a strawberry using household chemicals. Our hypothesis was right.
6. REPORT YOUR RESULTS. Using a lysis solution we can open the cells (that contain DNA) and thanks to the alcohol we produce precipitation, so the DNA that is in solution becomes solid again. Therefore, we can see and touch it.

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